

Product Information

DLPC

1,3-bis(sn-3'-phosphatidyl)-sn-glycerol,LP-R4-101

Cat. No.: **MPL0003K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

Phosphatidylcholine (PC)

Appearance

White powder

CAS Number

18194-25-7

Molecular Weight

621

Solubility

Chloroform

Purity

>99%

Size

500mg

Storage

2-8°C

Storage Shelf Time

3 years