

Product Information

DOPC

1,3-dioleoyl-sn-glycero-3-phosphocholine

Cat. No.: **MPL0005K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

Phosphatidylcholine (PC)

Appearance

White powder

CAS Number

4235-95-4

Molecular Weight

786

Solubility

Chloroform

Purity

>99%

Size

500mg

Storage

2-8°C

Storage Shelf Time

3 years