

Product Information

DOPE

1,2-dioleoyl-sn-glycero-3-phosphoethanolamine

Cat. No.: **MPL0015K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

Phosphadylethanolamine (PE)

Appearance

White solid

CAS Number

4004-05-1

Molecular Weight

744

Solubility

Chloroform

Purity

>99%

Size

500mg

Storage

-20°C

Storage Shelf Time

3 years