

Product Information

DPPE-PEG-Alkyne 16:0

Cat. No.: **MPL0029K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

Phosphadylethanolamine (PE)

Appearance

White powder/solid

Molecular Weight

PEG molecular weight available: 350, 550, 750, 1000, 2000, 3400, 5000, 10000, 20000, 30000, 40000

Solubility

Chloroform; methanol; DMSO

Purity

>95%

Size

100mg

Storage

-20°C

Storage Shelf Time

1 year