

Product Information

DSPG

1,2-distearoyl-sn-glycero-3-phospho-(1'-rac-glycerol) (sodium salt)

Cat. No.: **MPL0089K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

Phosphatidylglycerol (PG)

Appearance

White powder

CAS Number

200880-42-8

Molecular Weight

801

Solubility

CHCl₃/MeOH/H₂O (65:35:8)

Purity

>99%

Size

500mg

Storage

-20°C

Storage Shelf Time

3 years