

Product Information

DPPG

1,2-dipalmitoyl-sn-glycero-3-phospho-(1-rac-glycerol)

Cat. No.: **MPL0090K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

Phosphatidylglycerol (PG)

Appearance

White powder

CAS Number

200880-41-7

Molecular Weight

744

Solubility

CHCl₃/MeOH/H₂O (65:35:8)

Purity

>99%

Size

500mg

Storage

-20°C

Storage Shelf Time

3 years