

Product Information

DLPA

1,2- Dilauroyl -sn-glycero-3-phosphate

Cat. No.: **MPL0096K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

Phosphatidic Acid (PA)

Appearance

White powder

CAS Number

108321-06-8

Molecular Weight

558

Solubility

CHCl₃/MeOH/H₂O (65:35:8)

Purity

>99%

Storage

-20°C

Storage Shelf Time

3 years