

Product Information

DOPA

1,2-dioleoyl-sn-glycero-3-phosphate (sodium salt)

Cat. No.: **MPL0097K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

Phosphatidic Acid (PA)

Appearance

White powder

CAS Number

108392-02-5

Molecular Weight

722

Solubility

Chloroform

Purity

>99%

Size

500mg

Storage

-20°C

Storage Shelf Time

6 months