

Product Information

DPPS

1,2-Dipalmitoyl-sn-glycero-3-phospho-L-serine

Cat. No.: **MPL0098K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

Phosphatidylserine (PS)

Appearance

White powder

CAS Number

145849-32-7

Molecular Weight

757

Solubility

CHCl₃/MeOH/H₂O (65:35:8)

Purity

>99%

Size

500mg

Storage

-20°C

Storage Shelf Time

3 years