

Product Information

DSPE-Mpeg2000

1,2-distearoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-2000]

Cat. No.: **MPL0104K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

White powder

CAS Number

147867-65-0

Molecular Weight

≈2800

Solubility

H₂O/Chloroform

Purity

>99%

Size

500mg

Storage

-20°C

Storage Shelf Time

3 years