

# Product Information

## DSPE-PEG2000-cRGD

**cRGD is an active targeting polypeptide. Liposomes carrying cRGDs can target tumor target.**

Cat. No.: **MPL0105K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PEGylation, Lipids (PEG-PE)

### Appearance

White solid

### Molecular Weight

PEG:2000 Da

### Solubility

H<sub>2</sub>O/Dichloromethane

### Purity

>90%

### Size

100mg

### Storage

-20°C

### Storage Shelf Time

1 year