

Product Information

Cyanine3-PEG-DSPE,MW:2000

DSPE-PEG-Cyanine3

Cat. No.: **MPL0119K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

Red Powder

Molecular Weight

PEG2000

Solubility

H₂O

Purity

>90%

Size

100mg

Storage

-20°C, avoiding light

Storage Shelf Time

1 year