

Product Information

DSPE-PEG-Cyanine5,MW:2000

Cat. No.: **MPL0121K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

Pale blue powder or solid

Molecular Weight

3100

Solubility

Dichloromethane:Methanol = 1:1

Purity

>90%

Size

100mg

Storage

-20°C

Storage Shelf Time

1 year