

Product Information

DSPE-PEG-Cyanine5.5, MW:2000Da

Cat. No.: **MPL0122K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

Dark blue powder

Molecular Weight

PEG:2000

Solubility

H₂O/DMSO

Purity

>90%

Size

100mg

Storage

-20°C, avoiding light

Storage Shelf Time

1 year