

Product Information

DSPE-PEG-Rhodamine B, MW:2000

DSPE-PEG2000 labeled red dye Rhodamine B can be used as fluorescent liposomes.

Cat. No.: **MPL0128K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

Red solid

Solubility

Chloroform

Purity

>90%

Size

100mg

Storage

-20°C

Storage Shelf Time

1 year