

Product Information

DSPE-PEG-cRGD,MW:2000

Cat. No.: **MPL0144K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

White solid

Molecular Weight

PEG:550,750,1000,2000,3400,5000Da

Solubility

Dichloromethane

Purity

>90%

Size

100mg

Storage

-20°C

Storage Shelf Time

1 year