

---

# Product Information

## DSPE-PEG-Metronidazole, MW:2000

Cat. No.: **MPL0167K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PEGylation, Lipids (PEG-PE)

### Appearance

White powder

### Molecular Weight

PEG2000

### Solubility

Dichloromethane

### Purity

>95%

### Size

100mg

### Storage

-20°C

### Storage Shelf Time

1 year