

# Product Information

## DSPE-PEG-Cyanine7,MW:2000

Cat. No.: **MPL0181K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PEGylation, Lipids (PEG-PE)

### Appearance

Yellow powder/Yellow liquid

### Molecular Weight

3466

### Solubility

Organic solvent

### Purity

>90%

### Size

100mg

### Storage

-20°C,avoiding light

### Storage Shelf Time

1 year