

# Product Information

## DSPE-PEG2K-WSW (WSWGPYSC)

Cat. No.: **MPL0191K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PEGylation, Lipids (PEG-PE)

### Appearance

White powder/solid

### Molecular Weight

PEG molecular weight available: 350, 550, 750, 1000, 2000, 3400, 5000, 10000, 20000, 30000, 40000

### Solubility

Chloroform/Methanol/DMSO

### Purity

>95%

### Size

100mg

### Storage

-20°C

### Storage Shelf Time

1 year