

Product Information

DSPE-PEG-Cys-Pro-Leu-Gly-Val-Arg-Lys-Biotin

Cat. No.: **MPL0257K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

White powder/solid

Molecular Weight

PEG molecular weight available: 350, 550, 750, 1000, 2000, 3400, 5000, 10000, 20000, 30000, 40000

Solubility

Chloroform/Methanol/DMSO

Purity

>95%

Size

100mg

Storage

-20°C

Storage Shelf Time

1 year