

# Product Information

## DSPE-TK-PEOz MW:2000

Cat. No.: **MPL0268K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PEGylation, Lipids (PEG-PE)

### Appearance

Powder or solid

### Molecular Weight

2000

### Solubility

Water; DMSO; DCM; DMF

### Purity

>95%

### Storage

-20°C, avoiding light

### Storage Shelf Time

1 year