

Product Information

DSPE-PEG-FA

Cat. No.: **MPL0279K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

Light yellow powder

CAS Number

1236288-25-7

Molecular Weight

200, 400, 600, 800, 1K, 2K, 3.4K, 5K, 10K, 20K

Solubility

Organic solvents

Purity

>95%

Storage

-20°C

Storage Shelf Time

1 year