

# Product Information

## **DSPE-PEG3400-c(RGDfk)**

Cat. No.: **MPL0560K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

### **Product Type**

Phospholipids

### **Classification**

PEGylation, Lipids (PEG-PE)

### **Molecular Weight**

PEG molecular weight available: 350, 550, 750, 1000, 2000, 3400, 5000, 10000, 20000, 30000, 40000

### **Solubility**

Organic solvents

### **Purity**

>95%

### **Storage**

-20°C

### **Storage Shelf Time**

1 year