

# Product Information

## HBVpreS1, 2-21myr-PEG2000-DSPE

Cat. No.: **MPL0640K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PEGylation, Lipids (PEG-PE)

### Molecular Weight

PEG molecular weight available: 350, 550, 750, 1000, 2000, 3400, 5000, 10000, 20000, 30000, 40000

### Solubility

Organic solvents

### Purity

>95%

### Storage

-20°C

### Storage Shelf Time

1 year