

Product Information

DPPE-PEG-Multi-Arms (4 arm 8 arm)

Cat. No.: **MPL0677K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

Powder or Solid

Molecular Weight

PEG molecular weight available: 350, 550, 750, 1000, 2000, 3400, 5000, 10000, 20000, 30000, 40000

Purity

>95%

Storage

-20°C

Storage Shelf Time

1 year