

---

# Product Information

## DLPE-PEG-Acrylamide

Cat. No.: **MPL0733K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PEGylation, Lipids (PEG-PE)

### Appearance

Powder or Solid

### Molecular Weight

PEG molecular weight available: 350, 550, 750, 1000, 2000, 3400, 5000, 10000, 20000, 30000, 40000

### Purity

>95%

### Storage

-20°C

### Storage Shelf Time

1 year