

---

## Product Information

### DLPE-PEG-Protected Amines (FMOC, tBOC)

Cat. No.: **MPL0743K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

#### Product Type

Phospholipids

#### Classification

PEGylation, Lipids (PEG-PE)

#### Appearance

Powder or Solid

#### Molecular Weight

PEG molecular weight available: 350, 550, 750, 1000, 2000, 3400, 5000, 10000, 20000, 30000, 40000

#### Purity

>95%

#### Storage

-20°C

#### Storage Shelf Time

1 year