

Product Information

DSPE-PEG-SH

Cat. No.: **MPL0817K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

Gray/white solid or viscous liquid depending on molecular weight

Molecular Weight

PEG molecular weight available: 350, 550, 750, 1000, 2000, 3400, 5000, 10000, 20000, 30000, 40000

Solubility

Most organic solvents or H₂O

Purity

>95%

Size

1 g

Storage

-20°C

Storage Shelf Time

1 year