

Product Information

DSPE-PEG-Mal

Cat. No.: **MPL0818K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

White solid

CAS Number

474922-22-0

Molecular Weight

PEG molecular weight available: 350, 550, 750, 1000, 2000, 3400, 5000, 10000, 20000, 30000, 40000

Solubility

Chloroform and warm water

Purity

>95%

Size

1 g

Storage

-20°C

Storage Shelf Time

1 year