

# Product Information

## DSPE-PEG-COOH

### DSPE-PEG-COOH

Cat. No.: **MPL0820K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PEGylation, Lipids (PEG-PE)

### Appearance

White solid

### CAS Number

1403744-37-5

### Molecular Weight

PEG molecular weight available: 2000, 3400, 5000

### Solubility

Chloroform and warm water

### Purity

>90%

### Size

500 mg; 1 g

### Storage

-20°C

### Storage Shelf Time

1 year