

Product Information

DSPE-PEG-Biotin

DSPE-PEG-Biotin

Cat. No.: **MPL0821K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

White solid

CAS Number

385437-57-0

Molecular Weight

PEG molecular weight available: 2000, 3400, 5000

Solubility

Organic solvents

Purity

>90%

Size

500 mg; 1 g

Storage

-20°C

Storage Shelf Time

1 year