

Product Information

DSPE-PEG-FITC

DSPE-PEG-FITC

Cat. No.: **MPL0823K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Molecular Weight

PEG molecular weight available: 2000, 3400, 5000

Solubility

Water; Methanol; Chloroform; DMF; DMSO

Purity

>90%

Size

100 mg; 500 mg

Storage

-20°C

Storage Shelf Time

1 year