

Product Information

DSPE-PEG-Galactose

DSPE-PEG-Galactose

Cat. No.: **MPL0824K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

White solid

Molecular Weight

PEG molecular weight available: 2000, 3400, 5000

Solubility

Methanol:Chloroform=1:1

Purity

>90%

Size

100 mg; 500 mg

Storage

-20°C

Storage Shelf Time

1 year