

Product Information

Cholesterol-PEG-OH

Cholesterol-PEG-OH

Cat. No.: **MPL0834K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

CLS-PEG

Molecular Weight

PEG molecular weight available: 2000, 3400, 5000

Solubility

DMSO, DMF, DCM

Purity

>95%

Size

500 mg; 1 g

Storage

-20°C

Storage Shelf Time

1 year