

# Product Information

## Cholesterol-PEG-Rhodamine

### CLS-PEG-RB; Cholesterol-PEG-Rhodamine

Cat. No.: **MPL0835K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

### Product Type

Phospholipids

### Classification

CLS-PEG

### Molecular Weight

PEG molecular weight available: 2000, 3400, 5000

### Solubility

Water; Methanol; Chloroform; DMSO

### Purity

>95%

### Size

100 mg; 500 mg

### Storage

-20°C

### Storage Shelf Time

1 year