

Product Information

Cholesterol-PEG-FA

Cholesterol-PEG-FA; CLS-PEG-FA

Cat. No.: **MPL0836K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

CLS-PEG

Appearance

Solid

Molecular Weight

PEG molecular weight available: 2000, 3400, 5000

Solubility

DMSO, DMF

Purity

>95%

Size

100 mg; 500 mg

Storage

-20°C

Storage Shelf Time

1 year