

---

# Product Information

## DPhPC

### 1,2-diphytanoyl-sn-glycero-3-phosphocholine

Cat. No.: **MPL0843K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

### Product Type

Phospholipids

### Classification

Phosphatidylcholine (PC)

### Appearance

White solid

### CAS Number

207131-40-6

### Molecular Weight

846.252

### Purity

>95%

### Storage

-20°C

### Storage Shelf Time

1 year