
Product Information

Fmoc-NH-PEG1-Propanoic acid

Fmoc-6-amino-4-oxahexanoic acid

Cat. No.: **MPL0867K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

White solid

CAS Number

1654740-73-4

Molecular Weight

355.38

Purity

>98%

Storage

2-8°C

Storage Shelf Time

1 year