

---

# Product Information

## **FITC-PEG-FA**

### **Fluorescein (polyethylene glycol) folic acid**

Cat. No.: **MPL0868K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

### **Product Type**

Phospholipids

### **Classification**

PEGylation, Lipids (PEG-PE)

### **Appearance**

Solid

### **Molecular Weight**

PEG molecular weight available: 1000, 2000, 5000, 6000, 8000

### **Solubility**

Water; DMSO

### **Purity**

>95%

### **Storage**

-20°C

### **Storage Shelf Time**

1 year