

Product Information

Boc-NH-PEG-NHS

BOC (polyethylene glycol) Succinimidyl Carboxyl Methyl ester

Cat. No.: **MPL0870K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

Solid or Semi-solid

Molecular Weight

PEG molecular weight available: 1000, 2000, 5000, 10000, 20000

Solubility

Water; Chloroform; DMSO

Purity

>95%

Storage

-20°C

Storage Shelf Time

1 year