

---

# Product Information

## MPEG-FITC

### mPEG Fluorescein

Cat. No.: **MPL0874K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

### Product Type

Phospholipids

### Classification

PEGylation, Lipids (PEG-PE)

### Appearance

Solid or Semi-solid

### Molecular Weight

PEG molecular weight available: 1000, 2000, 3400, 5000, 10000

### Solubility

Water; Chloroform; DMSO

### Purity

>95%

### Storage

-20°C

### Storage Shelf Time

1 year