

Product Information

mPEG-NHS

Methoxy (polyethylene glycol) SuccinimidylCarboxyl Methyl ester

Cat. No.: **MPL0881K**

This product is for research use only and is not intended for diagnostic use.

This product is a phospholipid for research use only.

Product Type

Phospholipids

Classification

PEGylation, Lipids (PEG-PE)

Appearance

Solid or Semi-solid

CAS Number

92451-01-9

Molecular Weight

PEG molecular weight available: 1000, 2000, 3400, 5000, 10000

Solubility

Water; Chloroform; DMSO

Purity

>95%

Storage

-20°C

Storage Shelf Time

1 year