

Product Information

Clodronate liposome (Anionic)

Cat. No.: **MPL0912K**

This product is for research use only and is not intended for diagnostic use.

The phospholipid bilayer is composed of phosphatidylcholine, phosphatidylglycerol and cholesterol. The concentration of clodronate disodium is 20 mg/mL.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

CAS Number

22560-50-5

Purity

>95%

Size

2 mL; 10 mL

Storage

2-8°C

Storage Shelf Time

6 months