

# Product Information

## PEGylated HSPC, CHOL, mPEG2000-DSPE Liposomes

Cat. No.: **MPL0925K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of PEGylated HSPC, CHOL and mPEG2000-DSPE (50:45:5). The mean particle size is 100 nm and the lipid concentration is 60mM.

### Product Type

Phospholipids

### Classification

Liposome

### Appearance

Liquid

### Purity

>95%

### Size

2 mL; 10 mL

### Storage

2-8°C

### Storage Shelf Time

4 months