

# Product Information

## DC-CHOL, DOPE Liposomes

Cat. No.: **MPL0929K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of DC-CHOL and DOPE (30:70). The mean particle size is 90- 110 nm and the lipid concentration is 5mM.

### Product Type

Phospholipids

### Classification

Liposome

### Appearance

Liquid

### Purity

>95%

### Size

2 mL; 10 mL

### Storage

2-8°C

### Storage Shelf Time

4 months