

Product Information

DPPC, CHOL, mPEG-DSPE liposomes

Cat. No.: **MPL0936K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of DPPC, CHOL and mPEG2000-DSPE (50:45:5). The mean particle size is 100 nm and the lipid concentration is 60mM.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL; 10 mL

Storage

2-8°C

Storage Shelf Time

4 months