

# Product Information

## PEGylated Liposome with Doxorubicin HCl (2.0mg, mL)

Cat. No.: **MPL0937K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of HSPC, CHOL and mPEG2000-DSPE (56.2:38.5:5.3). The mean particle size is 85 nm and the lipid concentration is 20mM.

### Product Type

Phospholipids

### Classification

Liposome

### Appearance

Liquid

### Purity

>95%

### Size

2 mL; 10 mL

### Storage

2-8°C

### Storage Shelf Time

4 months