

Product Information

PEGylated Liposome with Doxorubicin HCl (4.0mg, mL)

Cat. No.: **MPL0938K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of HSPC, CHOL and mPEG2000-DSPE (56.2:38.5:5.3). The mean particle size is 85 nm and the lipid concentration is 40mM.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL; 10 mL

Storage

2-8°C

Storage Shelf Time

4 months