

Product Information

PEGylated Liposome with Topotecan HCl (2.0mg, mL)

Cat. No.: **MPL0939K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of HSPC, CHOL and mPEG2000-DSPE (50:45:5). The mean particle size is 100 nm and the lipid concentration is 20mM.

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL; 10 mL

Storage

2-8°C

Storage Shelf Time

4 months