

Product Information

Control of PEGylated Liposome with Doxorubicin (2.0mg, mL)

Cat. No.: **MPL0942K**

This product is for research use only and is not intended for diagnostic use.

The liposome is composed of Hydrogenated Soy PC, Cholesterol and DSPE-PEG(2000) (57:38:5). The mean particle size is 80-100 nm and the lipid concentration is 21.62mM. It should be used as a negative control for PEGylated Liposomal with Doxorubicin (2.0mg/mL).

Product Type

Phospholipids

Classification

Liposome

Appearance

Liquid

Purity

>95%

Size

2 mL; 5 mL

Storage

2-8°C

Storage Shelf Time

4 months